

COMPATIBILITY CHART

		Filter Media											Housing						
		Cellulose Acetate	Nitrocellulose	Polyethersulfone	Nylon 66	PTFE (unlaminated)	PTFE (laminated)	PVDF Philic	RC	Glass Fiber (binder)	Glass Fiber (no binder)	Silver	Polycarbonate	Polyester	Modified Acrylic	Polysulfone	Polystyrene	Polypropylene	
Chemical		ca	nc	pes	ny	ptu	ptl	pvdff	rc	gfb	gfn	ag	pc	pet	ac	ps	pst	pp	
ACIDS	Acetic Acid 5%	R	R	R	R	R	R	R	R	T	R	R	R	R	N	R	R	R	
	Acetic Acid 10%	N	N	R	L	R	R	R	R	T	R	R	R	R	N	R	R	R	
	Acetic Acid, Glacial	N	N	R	N	R	R	R	R	N	R	R	L	NR	N	R	R	L	
	Boric Acid	R	R	T	L	R	R	T	T	T	T	R	R	R	N	R	R	R	
	Hydrochloric, 6N	L	N	R	N	R	R	L	N	N	R	R	R	L	N	R	R	T	
	Hydrochloric, Conc.	N	N	R	N	R	R	R	N	N	R	R	R	N	N	R	R	T	
	Hydrofluoric, 10%	N	N	T	N	R	R	L	N	N	R	R	T	T	T	T	T	R	
	Hydrofluoric, 35%	N	N	T	N	R	R	T	R	N	N	R	T	T	T	T	T	T	
	Nitric Acid, 6N	L	R	N	N	R	L	T	N	N	L	N	R	R	N	N	L	T	
	Nitric Acid, Conc.	N	N	N	N	R	N	R	N	N	L	N	R	N	N	N	N	T	
	Sulfuric Acid, 6N	L	R	T	N	R	L	R	L	N	R	N	R	R	N	N	N	T	
	Sulfuric Acid, Conc.	N	N	N	N	R	N	T	N	N	R	N	N	N	N	N	N	T	
	ALCOHOLS	Amyl Alcohol	R	N	N	R	R	R	R	R	R	R	R	T	T	N	R	N	R
		Benzyl Alcohol	L	R	N	L	R	R	R	R	N	N	R	NR	NR	R	R	N	R
Butyl Alcohol		R	R	R	R	R	R	R	T	R	R	R	R	R	R	R	T	R	
Butyl Cellosolve		L	N	T	R	R	R	T	T	R	R	R	L	R	T	L	T	T	
Ethyl Alcohol <80%		R	R	R	R	R	R	R	T	R	R	R	R	R	L	R	L	T	
Ethyl Alcohol >80%		R	L	R	R	R	R	R	T	R	R	R	R	R	L	R	N	T	
Ethylene Glycol		R	L	R	R	R	R	R	R	R	R	R	R	R	T	R	T	R	
Glycerine (Glycerol)		R	R	R	R	R	R	R	R	R	R	R	R	R	T	R	T	R	
Isobutyl alcohol		R	R	T	R	R	R	R	T	N	N	R	R	R	R	R	R	T	
Isopropanol		R	L	R	R	R	R	R	R	R	R	R	R	R	T	R	T	T	
Methanol		R	N	R	T	R	R	R	R	R	R	R	R	T	R	R	R	T	
Methyl Cellosolve		L	L	T	R	R	R	R	T	R	R	R	N	R	T	R	T	T	
Propanol		R	R	T	R	R	R	R	R	R	R	R	R	R	T	R	T	R	
BASES		Ammonium Hydroxide, 6N	N	N	R	N	R	R	L	N	R	R	N	L	R	R	R	T	
	Potassium Hydroxide, 6N	N	N	T	R	R	R	R	L	N	T	R	N	N	T	R	T	T	
	Sodium Hydroxide, 6N	N	N	R	N	R	R	R	L	N	T	R	N	NR	T	T	T	T	
SOLVENTS	Acetone	N	N	N	R	R	R	N	R	R	R	R	L	R	N	N	N	R	
	Acetonitrile	N	N	R	T	R	R	R	R	T	R	T	NR	T	N	N	N	R	
	Amyl Acetate	L	N	L	R	R	R	R	N	R	R	R	R	R	N	N	N	L	
	Aniline	N	N	R	R	R	R	T	R	T	T	R	N	R	T	N	N	L	
	Benezene	L	R	R	T	R	L	R	R	N	R	R	NR	R	N	N	N	L	
	Bromoform	N	R	T	R	R	R	T	T	R	R	R	N	R	T	N	T	T	
	Butyl Acetate	L	N	L	R	R	R	T	R	N	R	R	R	R	N	N	N	L	
	Carbon Tetrachloride	L	R	R	R	R	L	R	R	N	N	R	NR	R	N	N	N	N	
	Cellosolve	R	N	T	R	R	R	T	R	R	R	R	R	R	N	N	T	T	
	Chloroform	N	R	N	NR	R	L	R	R	R	R	R	N	R	N	L	N	L	
	Cyclohexane	R	R	T	R	R	R	T	R	R	R	R	R	R	N	R	T	R	
	Cyclohexanone	N	N	N	T	R	R	N	R	R	R	R	L	T	N	N	N	R	
	Diethyl Acetamide	N	N	T	R	R	N	T	R	R	R	R	NR	NR	N	N	N	T	
	Dimethyl Formamide	N	N	N	R	R	R	N	L	N	R	R	NR	NR	N	N	N	R	
	Dimethyl Sulfoxide (DMSO)	N	N	N	R	R	R	N	R	N	R	T	N	R	N	N	N	T	
	Dioxane	N	N	L	R	R	R	R	R	R	R	R	N	R	N	N	N	R	
	Ethyl Ether	L	L	R	R	R	R	R	R	T	R	R	R	R	N	L	N	N	
	Ethylene Dichloride	L	L	T	R	R	R	T	R	R	R	R	N	R	T	N	T	T	
	Formaldehyde	L	N	R	R	R	R	R	T	R	R	R	R	R	N	R	N	R	
	Freon TF	R	R	R	R	R	R	R	T	R	R	R	R	R	L	R	N	T	
	Gasoline	R	R	T	R	R	R	R	R	R	R	R	R	R	N	R	N	N	
	Hexane	R	R	T	R	R	R	R	R	L	R	R	R	R	N	R	N	T	
	Isopropyl Acetate	N	N	T	R	R	R	N	R	N	R	R	R	R	N	N	N	R	
	Kerosene	R	R	T	R	R	R	R	R	R	R	R	R	R	N	N	N	T	
	Methyl Acetate	N	N	T	R	R	R	R	R	N	R	R	N	R	N	N	N	R	
	Methyl Ethyl Ketone (MEK)	N	N	N	R	R	R	NR	R	R	R	R	NR	R	N	N	N	T	
	Methyl Isobutyl Ketone	N	N	T	R	R	R	N	R	R	R	R	NR	T	N	N	N	T	
	Methylene Chloride	N	N	N	T	R	R	R	NR	R	R	R	N	NR	N	N	N	N	
	Nitrobenzene	N	N	N	T	R	R	R	NR	N	N	T	N	NR	N	N	N	R	
	Pentane	R	R	R	R	R	L	R	NR	R	R	R	R	R	N	R	N	T	
	Perchloroethylene	R	R	N	R	R	R	T	R	N	N	R	T	T	N	L	N	L	
	Pyridine	N	N	N	T	R	R	N	R	N	R	R	N	T	N	N	N	L	
	Tetrahydrofuran	N	N	N	T	L	L	N	R	T	L	R	N	T	N	N	N	L	
	Toluene	L	R	N	R	R	L	R	R	N	R	R	L	R	N	N	N	L	
Trichloroethane	L	N	L	T	R	R	T	NR	T	T	R	N	T	N	N	N	T		
Trichloroethylene	R	R	R	T	L	L	R	R	N	N	R	N	R	N	N	N	N		
Triethylamine	R	L	T	R	R	R	T	R	R	R	R	L	R	T	N	T	T		
Xylene	R	R	L	T	R	L	R	R	R	R	R	NR	NR	N	N	N	R		
MISCELLANEOUS	Cottonseed Oil	R	R	T	R	R	R	T	T	L	R	R	R	T	R	R	T	R	
	Hydrogen Peroxide (30%)	R	R	T	R	R	R	R	R	R	R	R	R	R	T	R	T	R	
	Kodak KMER FTFR	N	N	T	R	R	R	T	T	N	N	R	R	R	N	R	N	T	
	Peanut Oil	R	R	T	R	R	R	T	T	R	R	R	R	R	T	R	T	T	
	Petroleum Oils	T	R	L	T	R	T	R	R	T	T	R	R	R	T	T	T	R	
	Sesame Oil	R	R	T	R	R	R	T	T	R	R	R	R	R	T	R	T	T	
	Shipley (AS-111,340,1350)	N	N	T	R	R	R	T	T	N	N	R	R	R	N	R	N	T	
	Silicone Oils	R	R	R	R	R	R	R	R	R	R	R	R	R	T	R	T	R	
	Turpentine	R	R	T	R	R	R	T	T	R	R	R	R	R	T	R	T	T	
Waycoat 59	N	N	T	R	R	R	T	T	N	N	R	R	R	N	R	N	T		



THE ONLY WAY TO SAY FILTRATION

Our technology supports life and guarantees safety in the most critical environments



GVS Filtration Sdn. Bhd.
 Lot No 10F-2B, 10th Floor, Tower 5, PFCC, Jalan Puteri 1/2,
 Bandar Puteri, 47100 Puchong Selangor, Malaysia.
 +603 7800 0062 | lifesciences.my@gvs.com

www.gvsmalaysia.com.my

Need more assistance?
 Simply scan and connect with us instantly.



To explore our full range

