




## 1+PAC FILTRATION MEMBRANE

Part number: 1216721

Filter, Disk, NC, 0.45µm  
White. Individually packed.  
Without Pad Gridded. Sterile  
47mm. S1. Pack size: 100/pk

GVS Mixed Cellulose Esters (MCE) Filtration Membrane is an unsupported, hydrophilic membrane. Its **rapid flow rate** and **high throughput** make it ideal for use in bioburden and sterility testing.

- Recovery and retention of E.coli in water/waste-water analysis
- Sterility & Bioburden testing in QC/QA laboratories



**SAVE 40%  
WHEN YOU BUY NOW**



\*Promotion is now applicable to purchases of ANY **1+PAC FILTRATION MEMBRANE**, part number: **1216721** indicate the "PAC40" promo code in your Purchase Order to enjoy the discounts instantly.

Promotion period : **1st Sept. 2022 - 30th Nov. 2022**